

Anti-Zbtb7A/Fbi-1/Lrf Rabbit Monoclonal Antibody For Wb, IHC-P, IF/ICC, Elisa - O95365

Item Number: CM0003390



Introduction

High-quality rabbit monoclonal antibody against human ZBTB7A (also known as Pokemon, LRF). Validated for Western blot, IHC-P, IF/ICC, ELISA. Ideal for studying transcriptional regulation and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	ZBTB7A/FBI-1/LRF Rabbit mAb
Remarks/Alias	LRF; FBI1; FBI-1; TIP21; ZBTB7; MNDLFH; ZNF857A; pokemon; ZBTB7A/FBI-1/LRF
Species	Human
Gene ID (Human)	51341
Gene ID	51341
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 450-550 of human ZBTB7A/FBI-1/LRF (O95365).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O95365
Protein Weight	61kDa
Shipping	Ice bag