

Anti-Zbtb8A Rabbit Polyclonal Antibody For Wb, Elisa - Q96Br9

Item Number: CM0023047

Introduction

Rabbit polyclonal antibody against ZBTB8A (BOZF1), validated for WB and ELISA in human, mouse, rat. Targets a zinc finger transcription factor involved in transcriptional regulation. Ideal for IVD assay development and research.

[Learn More](#)

Parameter	Value
Product Name	ZBTB8A Rabbit pAb
Remarks/Alias	BOZF1; ZBTB8; ZNF916A; ZBTB8A
Species	Human
Gene ID (Human)	653121
Gene ID	653121
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96BR9
Protein Weight	50kDa
Shipping	Ice bag