

Anti-Zc3H12A Polyclonal Antibody For Wb, Elisa - Q5D1E8

Item Number: CM0021359



Introduction

Rabbit polyclonal antibody against human ZC3H12A (MCPIP-1/Regnase-1), an endoribonuclease that modulates mRNA stability in immune responses. Validated for WB and ELISA, suitable for inflammation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ZC3H12A Rabbit pAb
Remarks/Alias	ZC3H12A; MCPIP; MCPIP-1; MCPIP1; Reg1; dj423B22.1; endoribonuclease ZC3H12A
Species	Human
Gene ID (Human)	80149
Gene ID	80149
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 339-599 of human ZC3H12A(NP_079355.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5D1E8
Protein Weight	65kDa
Shipping	Ice bag