

Anti-Zc3H12C Rabbit Polyclonal Antibody (Mcpip3) For Wb, If/Icc, Elisa - Q9C0D7

Item Number: CM0021641

Introduction

High-quality rabbit polyclonal antibody targeting human ZC3H12C (MCPIP3), validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for RNA regulation and endonuclease research.

[Learn More](#)

Parameter	Value
Product Name	ZC3H12C Rabbit pAb
Remarks/Alias	MCPIP3; ZC3H12C
Species	Human
Gene ID (Human)	85463
Gene ID	85463
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 735-883 of human ZC3H12C (NP_203748.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9C0D7
Protein Weight	99kDa
Shipping	Ice bag