

Anti-Zc3H12D Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - A2A288

Item Number: CM0021380

Introduction

ZC3H12D polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets probable ribonuclease and tumor suppressor, ~58 kDa. Suitable for cancer and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	ZC3H12D Rabbit pAb
Remarks/Alias	TFL; p34; MCPIP4; C6orf95; dj281H8.1; ZC3H12D
Species	Human
GeneID (Human)	340152
GeneID	340152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-230 of human ZC3H12D (NP_997243.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A2A288
Protein Weight	58kDa
Shipping	Ice bag