

Anti-Zc3H14 Polyclonal Antibody For Wb And Elisa - Q6Pjt7

Item Number: CM0017976

Introduction

Rabbit polyclonal antibody against ZC3H14, validated for WB and ELISA applications. Detects ZC3H14 in human, mouse, and rat samples. Suitable for circRNA biogenesis and RNA-binding protein research.

[Learn More](#)

Parameter	Value
Product Name	ZC3H14 Rabbit pAb
Remarks/Alias	SUT2, MRT56, UKp68, MSUT-2, NY-REN-37, ZC3H14
Species	Human
Gene ID	79882
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of human ZC3H14 (NP_997545.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6PJT7
Protein Weight	83kDa
Shipping	Ice bag