

Anti-Zc3H15 Polyclonal Antibody For Wb, Elisa - Q8Wu90

Item Number: CM0014205

Introduction

Rabbit polyclonal antibody targeting human ZC3H15 (DFRP1/LEREPO4), suitable for Western blot and ELISA. Cross-reacts with mouse. Useful for studying DRG1 regulation and GTPase activity.

[Learn More](#)

Parameter	Value
Product Name	ZC3H15 Rabbit pAb
Remarks/Alias	DFRP1; HT010; LEREPO4; MSTP012; ZC3H15
Species	Human
GeneID (human)	55854
GeneID	55854
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ZC3H15 (NP_060941.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WU90
Protein Weight	49kDa
Shipping	Ice bag