


Anti-Zc3H7A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iwr0

Item Number: CM0017802



Introduction

Anti-ZC3H7A rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat ZC3H7A, a miRNA biogenesis regulator. Ideal for nuclear protein research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	ZC3H7A Rabbit pAb
Remarks/Alias	ZC3H7; HSPC055; ZC3HDC7; ZC3H7A
Species	Human
Gene ID	29066
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ZC3H7A (NP_054872.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IWR0
Protein Weight	111kDa
Shipping	Ice bag