

Anti-Zcchc10 Rabbit Polyclonal Antibody For Elisa - Q8Tbk6

Item Number: CM0020327

Introduction

Anti-ZCCHC10 rabbit polyclonal antibody validated for ELISA with mouse cross-reactivity. Recombinant immunogen. ZCCHC10 is a zinc finger CCHC domain-containing protein involved in metal binding. Suitable for research applications.

[Learn More](#)

Parameter	Value
Product Name	ZCCHC10 Rabbit pAb
Remarks/Alias	ZCCHC10
Species	Human
Gene ID	54819
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human ZCCHC10 (NP_001287745.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	ELISA
Cross Reactivity	Mouse
SWISS	Q8TBK6
Protein Weight	21kDa
Shipping	Ice bag