

Anti-Zcchc11 Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Tax3

Item Number: CM0030247



Introduction

Rabbit polyclonal antibody targeting ZCCHC11 (TUT4) for Western blot, IF/ICC, and ELISA. Validated in human, mouse, and rat samples. Ideal for mRNA uridylation and miRNA biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ZCCHC11 Rabbit pAb
Remarks/Alias	PAPD3; TENT3A; ZCCHC11
Species	Human
GeneID (Human)	23318
GeneID	23318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human ZCCHC11 (NP_001009881.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5TAX3
Protein Weight	185kDa
Shipping	Ice bag