

Anti-Zcchc3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nud5

Item Number: CM0020382

Introduction

Rabbit polyclonal antibody targeting human ZCCHC3, validated in Human and Mouse for WB, IF/ICC, and ELISA. ZCCHC3 is a nucleic acid-binding protein involved in innate antiviral immunity.

[Learn More](#)

Parameter	Value
Product Name	ZCCHC3 Rabbit pAb
Remarks/Alias	C20orf99; ZCCHC3
Species	Human
Gene ID (Human)	85364
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NUD5
Protein Weight	44kDa
Shipping	Ice bag