

Anti-Zcchc7 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N3Z6

Item Number: CM0024067

Introduction

ZCCHC7 rabbit polyclonal antibody validated for WB and ELISA in human, mouse, and rat samples. Detects 63 kDa zinc-finger protein localized to nucleolus; recombinant immunogen for specific detection.

[Learn More](#)

Parameter	Value
Product Name	ZCCHC7 Rabbit pAb
Remarks/Alias	AIR1; HSPC086; ZCCHC7
Species	Human
GeneID (Human)	84186
GeneID	84186
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-231 of human ZCCHC7 (NP_115602.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8N3Z6
Protein Weight	63kDa
Shipping	Ice bag