


Anti-Zcwpw1 Polyclonal Antibody For Wb, Elisa - Q9H0M4

Item Number: CM0024849



Introduction

Rabbit polyclonal antibody targeting human ZCWPW1, validated for WB and ELISA, with mouse cross-reactivity. Suitable for investigating PRDM9-mediated histone methylation, meiotic DSB repair, and male fertility. Predicted molecular weight 72 kDa.

[Learn More](#)

Parameter	Value
Product Name	ZCWPW1 Rabbit pAb
Remarks/Alias	ZCW1; ZCWPW1
Species	Human
Gene ID	55063
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 389-648 of human ZCWPW1 (NP_060454.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H0M4
Protein Weight	72kDa
Shipping	Ice bag