

Anti-Zdhhc15 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Mv8

Item Number: CM0019924

Introduction

Rabbit polyclonal antibody targeting human palmitoyltransferase ZDHHC15, suitable for Western blotting and ELISA. Useful for investigating protein palmitoylation, endosomal trafficking, and synaptic function.

[Learn More](#)

Parameter	Value
Product Name	ZDHHC15 Rabbit pAb
Remarks/Alias	ZDHHC15; DHHC15; MRX91; palmitoyltransferase ZDHHC15
Species	Human
Gene ID	158866
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 231-331 of human ZDHHC15 (NP_659406.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96MV8
Protein Weight	17kDa/38kDa/39kDa
Shipping	Ice bag