

Anti-Zdhhc18 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nue0

Item Number: CM0017745

Introduction

Rabbit polyclonal antibody targeting ZDHHC18 (DHHC18), suitable for Western blot, immunohistochemistry (paraffin), and ELISA. Detects human, mouse, rat; ideal for studying palmitoylation and innate immunity.

[Learn More](#)

Parameter	Specification
Product Name	ZDHHC18 Rabbit pAb
Remarks/Alias	DHHC18; DHHC-18; ZDHHC18
Species	Human
GeneID (Human)	84243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-388 of human ZDHHC18 (NP_115659.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NUE0
Protein Weight	42kDa
Shipping	Ice bag