

Anti-Zdhhc2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9UIj5

Item Number: CM0020969

Introduction

High-quality anti-ZDHHC2 rabbit polyclonal antibody validated for WB and ELISA applications. Detects human and rat ZDHHC2, a palmitoyltransferase implicated in cancer metastasis and synaptic plasticity. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	ZDHHC2 Rabbit pAb
Remarks/Alias	DHHC2; ZNF372; ZDHHC2
Species	Human
Gene ID (Human)	51201
Gene ID	51201
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 270-367 of human ZDHHC2 (NP_057437.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UIJ5
Protein Weight	42kDa
Shipping	Ice bag