

Anti-Zdhhc20 Rabbit Polyclonal Antibody For Wb And Elisa - Q5W0Z9

Item Number: CM0024046



Introduction

Rabbit polyclonal antibody against human ZDHHC20, validated for WB and ELISA. This palmitoyltransferase antibody reacts with human, mouse, and rat samples, suitable for studying EGFR signaling, SARS-CoV-2 spike lipidation, and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	ZDHHC20 Rabbit pAb
Remarks/Alias	DHHC20; DHHC-20; 4933421L13Rik; ZDHHC20
Species	Human
GeneID (human)	253832
GeneID	253832
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 270-365 of human ZDHHC20 (NP_001316988.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5W0Z9
Protein Weight	42kDa
Shipping	Ice bag