

Anti-Zdhhc4 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Npg8

Item Number: CM0021650

Introduction

Rabbit polyclonal antibody targeting human ZDHHC4 (palmitoyltransferase), suitable for Western blot and ELISA. Detects human, mouse, and rat ZDHHC4 with a predicted molecular weight of 40 kDa. Ideal for palmitoylation and signaling research.

[Learn More](#)

Parameter	Value
Product Name	ZDHHC4 Rabbit pAb
Remarks/Alias	ZNF374; ZDHHC4
Species	Human
Gene ID (Human)	55146
Gene ID	55146
Immunogen	Synthetic peptide (amino acids 110-210 of human ZDHHC4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPG8
Protein Weight	40kDa
Shipping	Ice bag