

Anti-Zdhhc5 Rabbit Pab For Wb, Elisa - Q9C0B5

Item Number: CM0019477

Introduction

Rabbit polyclonal antibody against ZDHHC5, validated for WB and ELISA. Detects human and mouse ZDHHC5 (78 kDa), a palmitoyltransferase regulating cell adhesion, immunity, and fatty acid uptake. Ideal for IVD research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	ZDHHC5 Rabbit pAb
Remarks/Alias	DHHC5, ZNF375
Species	Human
GeneID	25921
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 130-230 of human ZDHHC5 (NP_056272.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9C0B5
Protein Weight	78kDa
Shipping	Ice bag