

Anti-Zdhhc8 Rabbit Polyclonal Antibody For Wb, Elisa - Q9ULc8

Item Number: CM0023043

Introduction

Rabbit polyclonal antibody to ZDHHC8, a palmitoyltransferase. Validated for WB and ELISA applications. Reacts with human and mouse. Ideal for research on protein palmitoylation and cholesterol efflux.

[Learn More](#)

Parameter	Value
Product Name	ZDHHC8 Rabbit pAb
Remarks/Alias	DHHC8; ZNF378; ZHHCL1; ZDHHC8
Species	Human
Gene ID (Human)	29801
Gene ID	29801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-380 of human ZDHHC8 (NP_037505.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULC8
Protein Weight	81kDa
Shipping	Ice bag