


Anti-Zeb1 Monoclonal Antibody For Wb, Elisa - P37275

Item Number: CM0011588



Introduction

Rabbit monoclonal antibody against human ZEB1, validated for WB and ELISA. Suitable for EMT and transcription regulation studies, targeting the zinc finger E-box-binding homeobox 1 repressor.

[Learn More](#)

Parameter	Value
Product Name	ZEB1 Rabbit mAb
Remarks/Alias	BZP; TCF8; AREB6; FECD6; NIL2A; PPCD3; ZFHEP; ZFH1A; DELTAEF1
Species	Human
Gene ID (Human)	6935
Gene ID	6935
Immunogen	Recombinant protein Recombinant protein of human ZEB1.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P37275
Protein Weight	116 kDa/122 kDa/124 kDa
Shipping	Ice bag