


Anti-Zeb1 Monoclonal Antibody For Wb, Elisa - P37275

Item Number: CM0008574



Introduction

Rabbit monoclonal ZEB1 antibody validated for Western blot and ELISA. Targets transcription factor 8, a repressor of E-cadherin involved in epithelial-mesenchymal transition. Cross-reacts with human and rat ZEB1.

[Learn More](#)

Parameter	Specification
Product Name	ZEB1 Rabbit mAb
Remarks/Alias	BZP; TCF8; AREB6; FECD6; NIL2A; PPCD3; ZFHEP; ZFHX1A; DELTAEF1; ZEB1
Species	Human
GeneID	6935
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 800-900 of human ZEB1 (NP_110378.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P37275
Protein Weight	124kDa
Shipping	Ice bag