


Anti-Zeb1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P37275

Item Number: CM0024928



Introduction

Rabbit polyclonal antibody targeting human ZEB1 (Zinc finger E-box-binding homeobox 1), validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse, rat. Suitable for studying EMT, cancer progression, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	ZEB1 Rabbit pAb
Remarks/Alias	BZP; TCF8; AREB6; FECD6; NIL2A; PPCD3; ZFHEP; ZFH1A; DELTAEF1; ZEB1
Species	Human
GeneID (human)	6935
GeneID	6935
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37275
Protein Weight	124 kDa
Shipping	Ice bag