

Anti-Zeb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60315

Item Number: CM0029518



Introduction

Rabbit polyclonal antibody targeting ZEB2 (SIP1/ZFH1B), a transcriptional repressor of E-cadherin. Validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for EMT and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	ZEB2 Rabbit pAb
Remarks/Alias	SIP1; SIP-1; ZFH1B; HSPC082; SMADIP1; ZEB2
Species	Human
GeneID (human)	9839
GeneID	9839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1050-1214 of human ZEB2 (NP_055610.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60315
Protein Weight	136kDa
Shipping	Ice bag