

Anti-Zfand3 Polyclonal Antibody For Ihc-P And Elisa - Q9H8U3

Item Number: CM0028052



Introduction

Anti-ZFAND3 polyclonal antibody (rabbit origin) against recombinant human ZFAND3 (aa 1-227). Reacts with human, mouse, and rat. Suitable for IHC-P and ELISA. UniProt Q9H8U3.

[Learn More](#)

Parameter	Value
Product Name	ZFAND3 Rabbit pAb
Remarks/Alias	TEX27; ZFAND3
Species	Human
GeneID (Human)	60685
GeneID	60685
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human ZFAND3 (NP_068762.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8U3
Protein Weight	25kDa
Shipping	Ice bag