


Anti-Zfp36 Rabbit Polyclonal Antibody For Wb, Elisa - P26651

Item Number: CM0018967



Introduction

High-quality rabbit polyclonal antibody targeting human ZFP36 (tristetraprolin), validated for WB and ELISA. Suitable for studying mRNA decay, inflammation, and tumor suppression pathways.

[Learn More](#)

Parameter	Value
Product Name	ZFP36 Rabbit pAb
Remarks/Alias	TTP; G0S24; GOS24; TIS11; NUP475; zfp-36; RNF162A; ZFP36
Species	Human
Gene ID (Human)	7538
Gene ID	7538
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-326 of human ZFP36 (NP_003398.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P26651
Protein Weight	34kDa
Shipping	Ice bag