

Anti-Zfpm1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ix07

Item Number: CM0019158



Introduction

Rabbit polyclonal antibody against ZFPM1 (FOG1), validated in WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Essential for erythroid and megakaryocytic differentiation via GATA factor interaction. Ideal for hematopoiesis and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	ZFPM1 Rabbit pAb
Remarks/Alias	FOG; FOG1; ZNF408; ZNF89A; ZC2HC11A; ZFPM1
Species	Human
Gene ID (Human)	161882
Gene ID	161882
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ZFPM1 (NP_722520.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IX07
Protein Weight	105kDa
Shipping	Ice bag