

Anti-Zfpm2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Ww38

Item Number: CM0027250

Introduction

High-quality rabbit polyclonal antibody targeting human ZFPM2 (FOG2). Validated for WB and ELISA with cross-reactivity to mouse. Ideal for studying transcription regulation, heart morphogenesis, and gonadal differentiation.

[Learn More](#)

Parameter	Value
Product Name	ZFPM2 Rabbit pAb
Remarks/Alias	DIH3, FOG2, SRXY9, ZNF89B, hFOG-2, ZC2HC11B, ZFPM2
Species	Human
Gene ID (Human)	23414
Gene ID	23414
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 982-1151 of human ZFPM2 (NP_036214.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WW38
Protein Weight	128kDa
Shipping	Ice bag