

Anti-Zfx Polyclonal Antibody For Wb - P17010

Item Number: CM0020792

Introduction

Rabbit polyclonal antibody against human ZFX, validated for Western blot (WB) and ELISA. Suitable for human, mouse, and rat samples. Zinc finger X-chromosomal protein, a probable transcriptional activator.

[Learn More](#)

Parameter	Value
Product Name	ZFX Rabbit pAb
Remarks/Alias	ZNF926; ZFX
Species	Human
GeneID (Human)	7543
GeneID	7543
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 180-380 of human ZFX (NP_003401.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17010
Protein Weight	91kDa
Shipping	Ice bag