


Anti-Zfyve9 Polyclonal Antibody For Wb And Elisa - 095405

Item Number: CM0024975



Introduction

Rabbit polyclonal antibody against human ZFYVE9 (SARA), suitable for WB and ELISA. Targets the early endosomal protein involved in TGF- β signaling and SMAD localization. Ideal for research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	ZFYVE9 Rabbit pAb
Remarks/Alias	NSP; SARA; MADHIP; SMADIP; PPP1R173; ZFYVE9
Species	Human
GeneID (Human)	9372
GeneID	9372
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ZFYVE9 (NP_004790.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	095405
Protein Weight	156kDa
Shipping	Ice bag