

## Anti-Zgpat Polyclonal Antibody For Wb, Elisa - Q8N5A5

Item Number: CM0024457



### Introduction

Rabbit polyclonal antibody targeting human ZGPAT for WB and ELISA applications. Shows cross-reactivity with mouse. ZGPAT is a transcriptional repressor with tumor suppressor potential.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | ZGPAT Rabbit pAb                                                                                                   |
| Remarks/Alias    | ZIP; ZC3H9; GPATC6; GPATCH6; ZC3HDC9; KIAA1847; ZGPAT                                                              |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 84619                                                                                                              |
| GeneID           | 84619                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 213-272 of human ZGPAT (NP_115916.3) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Mouse                                                                                                              |
| SWISS            | Q8N5A5                                                                                                             |
| Protein Weight   | 57kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |