

Anti-Zhx2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6X8

Item Number: CM0023674



Introduction

High-quality rabbit polyclonal antibody targeting human ZHX2, a transcriptional repressor involved in retinal development and neural progenitor maintenance. Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Specification
Product Name	ZHX2 Rabbit pAb
Remarks/Alias	RAF; AFR1; ZHX2
Species	Human
GeneID (Human)	22882
GeneID	22882
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 488-837 of human ZHX2 (NP_055758.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6X8
Protein Weight	92kDa
Shipping	Ice bag