

Anti-Zic3 Polyclonal Antibody For Wb And Elisa - O60481

Item Number: CM0027689



Introduction

CamelBio offers a ZIC3 rabbit polyclonal antibody validated for Western blot and ELISA, with cross-reactivity to human, mouse, and rat samples. Ideal for studying developmental processes and left-right asymmetry.

[Learn More](#)

Parameter	Value
Product Name	ZIC3 Rabbit pAb
Remarks/Alias	HTX; HTX1; ZNF203; VACTERLX; ZIC3
Species	Human
GeneID	7547
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ZIC3 (NP_003404.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60481
Protein Weight	51kDa
Shipping	Ice bag