


Anti-Zinc-A2-Glycoprotein (Zag/Azgp1) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P25311

Item Number: CM0002854



Introduction

Rabbit monoclonal antibody targeting human Zinc- α 2-glycoprotein (ZAG/AZGP1), validated for WB, IHC-P, ELISA. Ideal for studying lipid degradation and cancer-associated fat loss. UniProt: P25311.

[Learn More](#)

Parameter	Value
Product Name	Zinc- α 2-glycoprotein (ZAG/AZGP1) Rabbit mAb
Remarks/Alias	ZAG; ZA2G; Zinc- α 2-glycoprotein (ZAG/AZGP1)
Species	Human
Gene ID (Human)	563
Gene ID	563
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Zinc- α 2-glycoprotein (ZAG/AZGP1) (P25311).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P25311
Protein Weight	34 kDa
Shipping	Ice bag