


Anti-Zinc-A2-Glycoprotein (Zag/Azgp1) Rabbit Polyclonal Antibody For Wb, Elisa - P25311

Item Number: CM0028633



Introduction

Zinc- α 2-glycoprotein (ZAG/AZGP1) rabbit polyclonal antibody validated for WB and ELISA. Detects human and mouse ZAG, a secreted lipid-mobilizing factor implicated in cancer cachexia. Ideal for metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Zinc- α 2-glycoprotein (ZAG/AZGP1) Rabbit pAb
Remarks/Alias	ZAG; ZA2G; Zinc- α 2-glycoprotein (ZAG/AZGP1)
Species	Human
Gene ID (Human)	563
Gene ID	563
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25311
Protein Weight	34kDa
Shipping	Ice bag