

Anti-Zip12 Rabbit Polyclonal Antibody For Wb, Elisa - Q504Y0

Item Number: CM0021817

Introduction

Rabbit polyclonal antibody targeting the human ZIP12 zinc transporter. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studies of zinc transport, neurulation, and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	Zip12 Rabbit pAb
Remarks/Alias	ZIP12; ZIP-12; LZT-Hs8; bA570F3.1; Zip12
Species	Human
Gene ID	221074
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC39A12 (NP_001138667.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q504Y0
Protein Weight	77kDa
Shipping	Ice bag