

Anti-Zip14 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15043

Item Number: CM0015545

Introduction

High-quality rabbit polyclonal antibody targeting human ZIP14, validated in WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for metal ion transport and homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	ZIP14 Rabbit pAb
Remarks/Alias	HCIN; NET34; ZIP14; cig19; HMNDYT2; LZT-Hs4
Species	Human
GeneID	23516
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 101-200 of human ZIP14 (NP_001121903.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15043
Protein Weight	54kDa
Shipping	Ice bag