

Anti-Zip8 Rabbit Polyclonal Antibody For Wb And Elisa - Q9C0K1

Item Number: CM0024295

Introduction

Rabbit polyclonal antibody targeting human ZIP8 (SLC39A8), validated for Western blot (WB) and ELISA, cross-reacts with mouse and rat. Suitable for zinc/manganese transport and immunity research.

[Learn More](#)

Parameter	Value
Product Name	ZIP8 Rabbit pAb
Remarks/Alias	ZIP8; CDG2N; PP3105; BIGM103; LZT-Hs6
Species	Human
GeneID (Human)	64116
GeneID	64116
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 201-300 of human ZIP8 (NP_071437.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9C0K1
Protein Weight	50kDa
Shipping	Ice bag