


Anti-Zmat2 Polyclonal Antibody For Wb, Elisa - Q96Nc0

Item Number: CM0013282



Introduction

Rabbit anti-human ZMAT2 polyclonal antibody validated for WB and ELISA; cross-reacts with mouse and rat. Targets ZMAT2, a spliceosome component, suitable for pre-mRNA splicing studies.

[Learn More](#)

Parameter	Value
Product Name	ZMAT2 Rabbit pAb
Remarks/Alias	Snu23; Ptg-12; ZMAT2
Species	Human
Gene ID (Human)	153527
Gene ID	153527
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human ZMAT2 (NP_653324.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NC0
Protein Weight	24kDa
Shipping	Ice bag