

Anti-Zmiz2 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nf64

Item Number: CM0021803

Introduction

Rabbit polyclonal antibody against ZMIZ2, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Useful for studying nuclear receptor coactivation and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	ZMIZ2 Rabbit pAb
Remarks/Alias	NET27; ZIMP7; hZIMP7; TRAFIP20; ZMIZ2
Species	Human
Gene ID (Human)	83637
Gene ID	83637
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human ZMIZ2 (NP_001287888)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NF64
Protein Weight	97kDa
Shipping	Ice bag