


Anti-Zmpste24 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O75844

Item Number: CM0023826



Introduction

Rabbit polyclonal antibody against ZMPSTE24, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying lamin A processing, nuclear envelope dynamics, and antiviral mechanisms.

[Learn More](#)

Parameter	Value
Product Name	ZMPSTE24 Rabbit pAb
Remarks/Alias	HGPS; PRO1; FACE1; RSDM1; STE24; FACE-1; Ste24p; ZMPSTE24
Species	Human
GeneID (Human)	10269
GeneID	10269
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ZMPSTE24 (NP_005848.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75844
Protein Weight	55kDa
Shipping	Ice bag