


Anti-Zmynd11 Polyclonal Antibody For Wb, Elisa - Q15326

Item Number: CM0030428



Introduction

Rabbit polyclonal antibody recognizing human ZMYND11 (BS69/BRAM1), suitable for WB and ELISA, with cross-reactivity to mouse. A key chromatin reader and tumor suppressor for epigenetic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ZMYND11 Rabbit pAb
Remarks/Alias	BS69; BRAM1; MRD30; ZMYND11
Species	Human
GeneID (Human)	10771
GeneID	10771
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 278-567 of human ZMYND11 (NP_997644.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15326
Protein Weight	71kDa
Shipping	Ice bag