


Anti-Zmynd8 Polyclonal Antibody For Wb, Elisa - Q9Ulu4

Item Number: CM0024760



Introduction

Anti-ZMYND8 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human ZMYND8, a chromatin reader involved in transcriptional regulation, DNA repair, and tumor suppression. Ideal for epigenetics and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ZMYND8 Rabbit pAb
Remarks/Alias	RACK7; PRKCBP1; PRO2893; ZMYND8
Species	Human
Gene ID (Human)	23613
Gene ID	23613
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 959-1188 of human ZMYND8 (NP_898868.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9ULU4
Protein Weight	132kDa
Shipping	Ice bag