


Anti-Zmynd8 Polyclonal Antibody (Ko Validated) For Wb, Ip, Elisa - Q9Ulu4

Item Number: CM0021618



Introduction

KO validated rabbit polyclonal antibody targeting human ZMYND8. Validated for WB, IP, ELISA with human and mouse reactivity. Useful for chromatin reader, transcription regulation, and DNA repair research.

[Learn More](#)

Specification	Details
Product Name	[KO Validated] ZMYND8 Rabbit pAb
Remarks/Alias	RACK7; PRKCBP1; PRO2893; D8
Species	Human
GeneID (Human)	23613
GeneID	23613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 959-1188 of human ZMYND8 (NP_898868.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULU4
Protein Weight	132kDa
Shipping	Ice bag