

Anti-Znf10 Polyclonal Antibody For Wb, Elisa - P21506

Item Number: CM0028336

Introduction

High-quality ZNF10 Rabbit pAb for Western blot and ELISA, detects human Zinc finger protein 10 (ZNF10/KOX1), a putative transcriptional regulator. Suitable for chromatin and gene regulation research. SWISS: P21506.

[Learn More](#)

Parameter	Value
Product Name	ZNF10 Rabbit pAb
Remarks/Alias	KOX1; ZNF10
Species	Human
GeneID (human)	7556
GeneID	7556
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 474-573 of human ZNF10 (NP_056209.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P21506
Protein Weight	66kDa
Shipping	Ice bag