

Anti-Znf10 Polyclonal Antibody For Wb, Elisa - P21506

Item Number: CM0019062

Introduction

ZNF10 Rabbit pAb is a polyclonal antibody validated for WB and ELISA applications in human samples. It targets the zinc finger protein ZNF10 (KOX1), a transcription regulator with a molecular weight of 66kDa.

[Learn More](#)

Parameter	Value
Product Name	ZNF10 Rabbit pAb
Remarks/Alias	KOX1; ZNF10
Species	Human
GeneID (Human)	7556
GeneID	7556
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ZNF10 (NP_056209.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P21506
Protein Weight	66kDa
Shipping	Ice bag