

Anti-Znf124 (Hzf-16) Rabbit Polyclonal Antibody For Wb, Elisa - Q15973

Item Number: CM0015554

Introduction

Rabbit polyclonal antibody against ZNF124 (HZF-16), suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	ZNF124 Rabbit pAb
Remarks/Alias	ZK7; HZF16; HZF-16; ZNF124
Species	Human
Gene ID (Human)	7678
Gene ID	7678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-200 of human ZNF124 (NP_001284497.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q15973
Protein Weight	40kDa
Shipping	Ice bag