

Anti-Znf133 Rabbit Polyclonal Antibody For Wb, Elisa - P52736

Item Number: CM0014310

Introduction

Rabbit polyclonal antibody targeting ZNF133 (zinc finger protein 133), validated for WB and ELISA in human, mouse, and rat. Suitable for transcription regulation research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	ZNF133 Rabbit pAb
Remarks/Alias	ZNF150; pHZ-13; pHZ-66; ZNF133
Species	Human
Gene ID (Human)	7692
Gene ID	7692
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ZNF133 (NP_003425.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52736
Protein Weight	73kDa
Shipping	Ice bag