

Anti-Znf134 Rabbit Polyclonal Antibody For Wb And Elisa - P52741

Item Number: CM0023788

Introduction

Anti-ZNF134 rabbit polyclonal antibody validated for Western blot and ELISA. Targets human zinc finger protein 134, a putative transcriptional regulator. Suitable for research and diagnostic applications. UniProt P52741.

[Learn More](#)

Parameter	Value
Product Name	ZNF134 Rabbit pAb
Remarks/Alias	pHZ-15; ZNF134
Species	Human
Gene ID (Human)	7693
Gene ID	7693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-170 of human ZNF134 (NP_003426.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P52741
Protein Weight	48kDa
Shipping	Ice bag