

Anti-Znf141 Polyclonal Antibody For Ihc-P, Elisa - Q15928

Item Number: CM0014237

Introduction

Rabbit polyclonal antibody targeting human ZNF141, suitable for IHC-P and ELISA applications. Cross-reacts with mouse and rat. Useful for transcriptional regulation and limb development research.

[Learn More](#)

Parameter	Value
Product Name	ZNF141 Rabbit pAb
Remarks/Alias	D4S90; pHZ-44; ZNF141
Species	Human
GeneID (Human)	7700
GeneID	7700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 405-474 of human ZNF141 (NP_003432.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15928
Protein Weight	55kDa
Shipping	Ice bag